



Telephone: (845) 365-8200 * Website: www.dynarex.com
Toll Free: (888) DYNAREX * Website: www.thcnet.com
Fax: (845) 365-8201

MATERIAL SAFETY DATA SHEET

Page 1 of 2

Dtd. 4/4/13

I. PRODUCT IDENTIFICATION

PRODUCT/TRADE NAME: ALCOHOL PREP/SWAB – Reorder #'s AM-20200, 1103, 1104, 1106, 1113, 1114, 1116
1203, 1204.

HAZARD RATING (NFPA)

HEALTH: 1

FLAMMABILITY: 3

REACTIVITY: 0

SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-8200 (M-F DAYTIME) AT OTHER TIMES, CONTACT
THE LOCAL POISON CONTROL CENTER

CHEMICAL NAME: ISOPROPYL ALCOHOL

II. HAZARDOUS INGREDIENTS PER 29 CFR 1910.1200

Hazardous Ingredients	%	ACGIH TLV	CAS NUMBER
ISOPROPYL ALCOHOL	70	400 ppm	67-63-0

III. PHYSICAL/CHEMICAL CHARACTERISTICS

COLOR/ODOR/APPEARANCE: TOWELETTE/PAD/SWAB SATURATED WITH CLEAR LIQUID WITH ALCOHOL
ODOR..

BOILING POINT: N/A

FLASH POINT: 78°F

VAPOR DENSITY N/A

EVAPORATION RATE: N/A

SOLUBILITY IN WATER: COMPLETE

SPECIFIC GRAVITY (H₂O = 1):N/A

IV. FIRE & EXPLOSION HAZARD DATA

FLASH POINT (Method Used): 78°F TAG CLOSED CUP : LEL: 2.0 UEL: 12.0

EXTINGUISHING MEDIA: DRY CHEMICAL OR ALCOHOL TYPE FOAM, CARBON DIOXIDE

SPECIAL FIRE FIGHTING PROCEDURES: HANDLE AS FLAMMABLE LIQUID. USE RESPIRATORY PROTECTION
FOR

FIRE FIGHTING PERSONNEL

UNUSUAL FIRE AND EXPLOSION HAZARDS: CLASS 3 FLAMMABILITY

V. REACTIVITY DATA

STABILITY: STABLE

CONDITIONS TO AVOID NONE

INCOMPATIBILITY: NONE

HAZARDOUS DECOMPOSITION OR BYPRODUCT: NONE

POLYMERIZATION: WILL NOT OCCUR

CONDITIONS TO AVOID: NONE

Executive Headquarters: Ten Glenshaw Street, Orangeburg, NY 10962



Telephone: (845) 365-8200 * **Website:** www.dynarex.com
Toll Free: (888) DYNAREX * **Website:** www.thcnet.com
Fax: (845) 365-8201

ALCOHOL PREP/SWAB

Page 2 of 2

VI. HEALTH HAZARD DATA

EFFECTS OF OVER EXPOSURE:

SKIN: IF RASH APPEARS, DISCONTINUE USE

EYES: WILL CAUSE EYE STING IF SPLASHED

INHALATION: NONE

INGESTION NOT A NORMAL ROUTE OF ENTRY

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: DISCONTINUE USE IF RASH OR IRRITATION OCCURS

EYE CONTACT: FLUSH WITH COLD WATER FOR 15 MINUTES

INHALATION: REMOVE TO FRESH AIR

INGESTION: INDUCE VOMITING. CALL PHYSICIAN

TARGET ORGANS: N/A

VII. SPILL AND DISPOSAL PROCEDURE

SPILL CONTROL: ELIMINATE ALL SOURCES OF IGNITION

WASTE DISPOSAL METHOD: FLUSH SPILLS WITH WATER. HANDLE AS FLAMMABLE LIQUID. FOLLOW LOCAL, STATE AND FEDERAL REGULATIONS.

HANDLING AND STORAGE: STORE AWAY FROM HEAT AND SOURCES OF IGNITION

VIII. CONTROL MEASURES/PROTECTION

RESPIRATION: USE IF PERMISSABLE EXPOSURE LEVEL IS EXCEEDED WHEN HANDLING BULK LIQUID

VENTILATION: RECOMMENDED

PROTECTIVE GLOVES: NONE REQUIRED

EYE PROTECTION: ONLY IF SPLASHING IS EXPECTED

HYGIENIC PRACTICES: GOOD HOUSEKEEPING PRACTICES SHOULD BE FOLLOWED.

OTHER: DON'T ALLOW LARGE QUANTITIES OF WASTE TO ACCUMULATE

IX. TRANSPORT/SHIPPING

DOT SHIPPING NAME: N/A

TECHNICAL SHIPPING NAME: N/A

DOT SHIPPING CLASSIFICATION: DOT-NOT REGULATED

DOT ID NO.: N/A

DOT LABEL REQUIREMENTS: N/A

UN/NA NUMBER REGULATIONS: N/A

REPORTABLE QUANTITY: N/A

X. DISCLAIMER

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND DYNAREX CORPORATION ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.