



Telephone: (845) 365-8200 \* Website: www.dynarex.com  
 Toll Free: (888) DYNAREX \* Website: www.thcnet.com  
 Fax: (845) 365-8201

Reorder # 2500 & 6500 Series

Page 1 of 2  
 Reviewed on 4/11/13

## MATERIAL SAFETY DATA SHEET

Product Name: *Nitrile Exam Gloves*  
 CAS#: Not Established  
 LAT/IACO Hazard Rating: Not Restricted  
 D.O.T. Rating: Not Restricted

.....  
***Physical and Chemical Properties***

Re-order # 2511-2512-2513-2514

Nitrile Latex (NBR)  
 Potassium Hydroxide  
 Teric 320  
 Sulphur  
 ZDEC (Zinc diethyldithiocarbamate)  
 ZDBC (Zinc dibutyldithiocarbamate)  
 Zinc Oxide  
 Titanium Dioxide  
 Wingstay-L  
 BLU-BS(Color Pigment)  
 Color: Blue  
 Odor: N/A



Telephone: (845) 365-8200 \* Website: www.dynarex.com  
Toll Free: (888) DYNAREX \* Website: www.thcnet.com  
Fax: (845) 365-8201

***Normal Precautions***

Ventilation	Not necessary under normal conditions
Equipment	Not necessary
Spills/Clean Up	Not necessary

.....

***Fire and Explosion***

Stability	Stable
Flash Point	N/A
Extinguishing Media	Water spray/carbon dioxide /form.
Conditions to Avoid	None
Incompatible Materials	None

.....

***Other Precautions***

Storage	Inner boxes on cases and are to be stored in original packing in a cool, under room temperature away turn dust, humidity and direct sunlight.
Shelf life	If properly stored ,the Nitrile gloves should have a shelf life of minimum 24 months.

---

Prepared By: Dynarex Corp.  
Dated: 2011