



Telephone: (845) 365-8200 * Website: www.dynarex.com
Toll Free: (888) DYNAREX * Website: www.thcnet.com
Fax: (845) 365-8201

Reorder # 2300 & 2400 Series

Page 1 of 3
Reviewed on 1/30/14

MATERIAL SAFETY DATA SHEET

FOR LATEX EXAMINATION GLOVES MADE OF NATURAL RUBBER POWDERED AND POWDER-FREE

IDENTIFICATION

Product Name: Latex Examination Gloves, Sterile and Non-sterile
Chemical Name: Sulfur cross-linked cis 1,4 Poly-isoprene
Formula: Composite
Molecular Weight: Mixture

PHYSICAL DATA

Specific Gravity: 0.94
Taste: Nil
Appearance: Creamy White (powdered)
Light Yellowish (powder-free)
Odor: Sweetish (powdered)
Faint Chlorine (powder-free)

INGREDIENTS

	Material	Hazard
ZDEC	Zinc diethyl dithiocarbamate	May cause type IV hypersensitivity
ZMBT	Zinc mercaptobenzothiazole	May cause type IV hypersensitivity
KOH	Potassium hydroxide	May cause type IV hypersensitivity
S	Sulphur	May cause type IV hypersensitivity
ZnO	Zinc Oxide	May cause type IV hypersensitivity
TiO ₂	Titanium dioxide	May cause type IV hypersensitivity
Wingstay L	Butylated reaction product of p-cresol and dicyclopentadiene	May cause type I hypersensitivity

FIRE AND EXPLOSION HAZARD

Non-combustible at ambient temperatures, combustible when ignited.
No explosion hazard.

HEALTH HAZARD INFORMATION

Powdered gloves:

Irritant contact dermatitis may result from prolonged contact with the product, the causative agents being glove powder and surfactants and inadequate rinsing after use. There are no evidence of any adverse effects from available information of sulfur, zinc oxide and titanium dioxide. However, presence of residual ZDEC, ZMBT and anti-oxidants on the product surface may cause delayed type IV hypersensitivity e.g. allergic contact dermatitis, chemical allergy. Reactions include itching, burning sensations, blistering, reddening and pain. In chronic cases, users may develop dry and thickened skin, cracks, peeling and crusting.



Telephone: (845) 365-8200 * **Website:** www.dynarex.com
Toll Free: (888) DYNAREX * **Website:** www.thcnet.com
Fax: (845) 365-8201

Natural rubber latex contain a small quantity of natural proteins which may cause an immediate type I hypersensitive reactions such as urticaria and anaphylaxis in a small percentage of the population. Type I hypersensitivity is life threatening and affected users must seek immediate emergency treatment and to substitute with alternative low protein or non-latex containing products. Signs of type I hypersensitivity include initial itching, hives, swelling, runny nose, wheezing and breathing difficulty. Other symptoms may involve nausea, abdominal cramps, increased heart palpitation, decreased blood pressure and anaphylactic shock.

Powder-free gloves:

Irritant contact dermatitis may result from prolonged contact with the product, the causative agents being residual surfactants and certain processing chemicals. There are no evidence of any adverse effects from available information of sulfur, zinc oxide and titanium dioxide. However, presence of residual ZDEC, ZMBT and anti-oxidants on the product surface may cause delayed type IV hypersensitivity e.g. allergic contact dermatitis and chemical allergy. Reactions include itching, burning sensations, blistering, reddening and pain. In chronic cases, users may develop dry and thickened skin, cracks, peeling and crusting.

Natural rubber latex contain a small quantity of natural proteins which may cause an immediate type I hypersensitive reactions such as urticaria and anaphylaxis in a small percentage of the population. However, the process of chlorination to make the gloves powder-free also reduces the extractable protein on the glove surface to below 50µg per gram glove. At such low protein level, sensitization is extremely rare. In any event, a user who cannot tolerate latex gloves is advised to consider alternative non-latex gloves. Type I hypersensitivity is often accompanied by symptoms such as itching, hives, swelling, runny nose, wheezing and breathing difficulty. Other symptoms may involve nausea, abdominal cramps, increased heart palpitation, decreased blood pressure and anaphylactic shock.

General Information:

The product is intended as a hand protective device only and not as a toy which can cause suffocation. Keep them away from young children or people incapable of using them properly.

Users are advised against touching sensitive skin areas such as the eye, nose, lips and mouth after donning to minimize skin irritations on certain individuals. Because of its protective function, the product is unventilated except at the cuff. Users should replace the product frequently, wash and clean their hands prior to donning the product again to reduce skin fatigue and roughening on the areas in frequent contact with the product.

EMERGENCY PROCEDURE

SWALLOWING:

Attempt removal of product to clear airway and contact medical service immediately



Telephone: (845) 365-8200 * **Website:** www.dynarex.com
Toll Free: (888) DYNAREX * **Website:** www.thcnet.com
Fax: (845) 365-8201

SKIN IRRITATION:

Wash with soap and water

EYE/NOSE CONTACT:

Wash with facial cleanser and water

STORAGE INFORMATION

The product must be stored in its original packing in a well-ventilated area away from dust, sunlight and excessive heat. It must be kept away from unauthorized personnel, young children and those unfamiliar with its applications.

DISPOSAL RECOMMENDATION

We recommend that all used products be disposed off in a hygienic manner, either through incineration or buried in authorized landfills.

We also recommend that the used product be disposed off separately from domestic and other nonclinical wastes as they may be carriers of infectious agents.

Prepared By: Dynarex Corp.
Dated: 2011