



Telephone: (845) 365-8200 \* Website: www.dynarex.com  
Toll Free: (888) DYNAREX \* Website: www.thcnet.com  
Fax: (845) 365-8201

Reviewed on 4/4/13

Page 1 of 4

## MATERIAL SAFETY DATA SHEET

### I. Product/Company Identification

PRODUCT/TRADE NAME: Oral Swabsticks— **Reorder # 1217-1218-1219-1220**

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-8200 (M-F DAYTIME) AT OTHER TIMES, CONTACT THE LOCAL POISON CONTROL CENTER

### II. Hazardous Ingredients/Identity Information

Hazardous Components:	Concentration	CAS NO.
None		

### III Hazards Identification

#### HARZARD CLASSIFICATION

Skin:	No adverse effects expected.
Ingestion:	May be harmful if swallowed in large quantities.
Inhalation:	No adverse effects expected.
Eyes:	May cause eye irritation on direct contact.

### IV. First Aid Measures

Skin:	Does not apply
Ingestion:	If large quantity is ingested, drink 1-2 glasses of water. Get medical attention. Do not give anything by mouth to an unconscious person.
Eyes:	Flush eyes under running water for 15minutes. Get medical attention if irritation persists.
Inhalation:	Does not apply



Telephone: (845) 365-8200 \* Website: [www.dynarex.com](http://www.dynarex.com)  
Toll Free: (888) DYNAREX \* Website: [www.thcnet.com](http://www.thcnet.com)  
Fax: (845) 365-8201

## v. Fire-fighting Measures

Fire hazard: Inflammable  
Explosion Limits: None  
Unusual fire and explosion hazards: Not considered to be an explosion hazard.  
Suitable fire extinguishing Media: Use any means suitable for extinguishing surrounding fire.  
Special Protection Equipment: In the event of fire, wear full protective clothing and self-contained breathing apparatus with full face piece operated in the positive pressure mode.

## vi. Accidental Release Measures

Personal Precautions: N/A  
Methods for Cleaning up: N/A

## vii. Handling and Storage

Handling: See personal protection, Section VIII  
Storage: Store at controlled room temperature, 0-25°C  
Containers: Keep in original container.  
Incompatibilities: None known  
Other Precautions: Prohibit open flame, spark and smoking. Do not contact with strong oxidizer.

## viii. Exposure Controls / Personal Protection.

General Protective Measures: Avoid eye contact  
General ventilation: Normal room ventilation recommended  
Eye Protection: Use of safety goggles is recommended.  
Hand Protection: Not necessary  
Skin Protection: Not necessary



Telephone: (845) 365-8200 \* Website: www.dynarex.com  
Toll Free: (888) DYNAREX \* Website: www.thcnet.com  
Fax: (845) 365-8201

## IX. Physical and Chemical Properties.

Appearance:	Homogenous Clear Liquid/Gel
Color:	Colorless to white
Odor:	Mint Odor
Vapor Pressure:	N/A
Boiling Point:	120□
PH:	7.5
Viscosity:	TBD
Specific gravity:	1.18

## X. Stability and Reactivity.

Stability:	Stable
Conditions to avoid:	High temperature
Material to avoid:	Avoid contact with strong oxidizing agents, acids, bases, amines.

Hazardous Decomposition Products: None known

Hazardous Polymerization Products: Hazardous polymerization will not occur.

## XI. Toxicological Information

This product has not been tested as a whole. The formula was reviewed taking into consideration available safety-related information on the individual ingredients, similar ingredients, similar formulas, and is determined to be safe for its intended use.

## XII. Ecological Information.

No specific ecological data are available for this product. Refer to Section VI for environmental precautions.

## XIII. Disposal Considerations

All disposal procedures must comply with local, state, and federal laws and regulations (contact local or state environmental agency for specific rules.)

Do not dump into sewers, grounds or body of water.

**Executive Headquarters:** Ten Glenshaw Street, Orangeburg, NY 10962



**Telephone:** (845) 365-8200 \* **Website:** [www.dynarex.com](http://www.dynarex.com)  
**Toll Free:** (888) DYNAREX \* **Website:** [www.thcnet.com](http://www.thcnet.com)  
**Fax:** (845) 365-8201

**xiv. Transport Information**

Not regulated as DOT hazardous material.

**xv. Regulatory Information**

CLASSIFICATION IN ACCORDANCE WITH EEC DIRECTIVES.

Hazard Classification Symbol: N/A  
Risk Statements: N/A  
Safety Statements: N/A

**xvi. DISCLAIMER**

**CAUTION: PRODUCT FOR PROFESIONAL USE ONLY.**

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND DYNAREX CORPORATION ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.

Dtd. 10/28/2010